## Ácido Úrico

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Methods List	EndPoint
<b>Test</b>	
Description	Uric Acid
Test ID	UA
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	505
Filter 2	-
<b>Volume (µL)</b>	
Sample	10
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	600
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	2.3
Normal high	8.2
Linearity limit	20.0
Absorbance limit	

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Methods List	EndPoint
<b>Test</b>	
Description	Albumin
Test ID	ALBUM
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	546
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	600
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	G/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	35
Normal high	50
Linearity limit	60
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	Kinetic
<b>Test</b>	
Description	a-Amylase
Test ID	AMYL
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	405
Filter 2	-
<b>Volume (µL)</b>	
Sample	10
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	60
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	5670*
Multiple	-
Specific	-
(Calibration)	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
<b>Options</b>	
Normal low	0
Normal high	100
Linearity limit	2000
Absorbance limit	

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Methods List	Differential
<b>Test</b>	
Description	Bilirubin Direct
Test ID	Bil D
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	546
Filter 2	-
<b>Volume (µL)</b>	
Sample	50
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	300
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	2
Minimum conc.	-
Replicate Blank	On
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	
<b>Options</b>	
Normal low	0.0
Normal high	0.3
Linearity limit	10.0
Absorbance limit	

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Methods List	Differential
<b>Test</b>	
Description	Bilirubin Total
Test ID	Bil T
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	546
Filter 2	-
<b>Volume (µL)</b>	
Sample	13
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	300
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	On
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	
<b>Options</b>	
Normal low	0.1
Normal high	1.2
Linearity limit	30.0
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Calcium AS
Test ID	Ca AS
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	630
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	8,6
Normal high	10.3
Linearity limit	25.0
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardenha, 110

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Methods List	Kinetic
<b>Test</b>	
Description	CK-MB
Test ID	CK MB
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	20
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	8254*
Multiple	-
Specific	-
<i>(Calibration)</i>	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	
<b>Options</b>	
Normal low	1
Normal high	24
Linearity limit	1000
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	Kinetic
<b>Test</b>	
Description	CK-NAC
Test ID	CKNAC
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	20
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	180
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	4127*
Multiple	-
Specific	-
<i>(Calibration)</i>	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
<b>Options</b>	
Normal low	24
Normal high	190
Linearity limit	1000
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Cholesterol
Test ID	CHOL
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	505
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	600
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	0
Normal high	200
Linearity limit	750
Absorbance limit	

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Methods List	FixedTime
<b>Test</b>	
Description	Creatinine
Test ID	CREAT
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	505
Filter 2	-
<b>Volume (µL)</b>	
Sample	25
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	60
2nd incubation	120
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	
Decr.	-
Incr.	On
Function	
<b>Options</b>	
Normal low	0.6
Normal high	1.3
Linearity limit	15.0
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	Differential
<b>Test</b>	
Description	Fe
Test ID	Fe
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	578
Filter 2	-
<b>Volume (µL)</b>	
Sample	50
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	600
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	µg/dl
Decimal point	0
Minimum conc.	-
Replicate Blank	On
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
<b>Options</b>	
Normal low	23
Normal high	168
Linearity limit	1000
Absorbance limit	

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Methods List	Kinetic
<b>Test</b>	
Description	Alk Phosph
Test ID	ALP
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	405
Filter 2	-
<b>Volume (µL)</b>	
Sample	10
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	60
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	2757*
Multiple	-
Specific	-
(Calibration)	-
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
	-
<b>Options</b>	
Normal low	42
Normal high	119
Linearity limit	700
Absorbance limit	Kinetic

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Glucose
Test ID	GLUC
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	505
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	600
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	70
Normal high	115
Linearity limit	400
Absorbance limit	

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Methods List	Differential
<b>Test</b>	
Description	Lactate
Test ID	Lac
Manufacturer	DiaSys
Expire	-
Note	Serum or CSF
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	7
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	300
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	On
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	
<b>Options</b>	
Normal low	4.5
Normal high	19.8
Linearity limit	120.0
Absorbance limit	

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Methods List	Kinetic
<b>Test</b>	
Description	LDH-P
Test ID	LDH
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	120
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	16030*
Multiple	-
Specific	-
<i>(Calibration)</i>	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	On
Incr.	-
<b>Options</b>	
Normal low	225
Normal high	450
Linearity limit	1200
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	FixedTime
<b>Test</b>	
Description	Lipase
Test ID	LIPAS
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	578
Filter 2	-
<b>Volume (µL)</b>	
Sample	10
Reagent 1	500
Reagent 2	125
<b>Reading parameter</b>	
1 st incubation	120
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
<b>Options</b>	
Normal low	0
Normal high	60
Linearity limit	300
Absorbance limit	-

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Magnesium
Test ID	Mg
Manufacturer	DiaSys
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	546
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	1
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	1.8
Normal high	2.6
Linearity limit	5.0
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Total Protein
Test ID	TProt
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	546
Filter 2	-
<b>Volume (µL)</b>	
Sample	10
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	300
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	g/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	66
Normal high	88
Linearity limit	150
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardenha, 110

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Methods List	Kinetic
<b>Test</b>	
Description	AST
Test ID	AST
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	50
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	120
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	1745*
Multiple	-
Specific	-
<i>(Calibration)</i>	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	On
Incr.	-
<b>Options</b>	
Normal low	0
Normal high	35
Linearity limit	260
Absorbance limit	

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Methods List	Kinetic
<b>Test</b>	
Description	ALT
Test ID	ALT
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	50
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	120
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	U/l
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	1745*
Multiple	-
Specific	-
<i>(Calibration)</i>	
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	On
Incr.	-
<b>Options</b>	
Normal low	0
Normal high	41
Linearity limit	260
Absorbance limit	

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Methods List	EndPoint
<b>Test</b>	
Description	Triglycerids
Test ID	TRIGL
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	505
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	-
<b>Reading parameter</b>	
1 st incubation	600
2nd incubation	-
Stability	1800
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
Function	-
<b>Options</b>	
Normal low	0
Normal high	200
Linearity limit	1000
Absorbance limit	

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Methods List	FixedTime
<b>Test</b>	
Description	UREA
Test ID	UREA
Manufacturer	DDS
Expire	-
Note	Serum or plasma
<b>Wavelength (nm)</b>	
Filter 1	340
Filter 2	-
<b>Volume (µL)</b>	
Sample	5
Reagent 1	500
Reagent 2	
<b>Reading parameter</b>	
1 st incubation	40
2nd incubation	60
Stability	-
Sample replicate	1
<b>Result</b>	
Measure units	mg/dl
Decimal point	0
Minimum conc.	-
Replicate Blank	-
<b>Calibration</b>	
Kfactor	-
Multiple	On
Specific	-
(Calibration)	Daily
N.standard	1
Replicate	1
Offset	-
Concentration 1	-
Decr.	-
Incr.	On
<b>Options</b>	
Normal low	17
Normal high	43
Linearity limit	300
Absorbance limit	-