

APPLICATION NOTE
COBAS 6000
α2 – Macroglobulin AUT KIT

1. Reagent preparation

Sample: Ready for use
 Reagent 1: Buffer, ready for use
 Reagent 3: Antiserum, ready for use
 Calibrator: Ready for use

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[700] [340]					
Sample Volume Norm.	[20] [10] [180]					
Sample Volume Dec.	[20] [10] [180]					
Sample Volume Inc.	[2] [0] [0]					
Configuration:						
R1	[150] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[15] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/dL]					
Decimal Places	[1]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[2]	[2]	[2]	[5]	[20]	[20]
Diluted S. Volume	[0]	[5]	[5]	[5]	[5]	[10]
Diluent Volume	[0]	[158]	[78]	[95]	[180]	[180]

* Entered by operator

3. Ordering information

AMG/AUT-000 1 x 10 mL Antiserum
 5 x 25 mL Buffer
 MPS/STH-001 Protein Standard High , 1 mL
 139F003 Immunology Control Low, 1mL
 139F002 Immunology Control High, 1mL

APPLICATION NOTE

COBAS 6000

α1 – Microglobulin AUS KIT

1. Reagent preparation

Sample: Centrifuged urine
 Reagent 1: Buffer, ready for use
 Reagent 3: Latex, ready for use
 Calibrator: AMI Standard, ready for use. Use 9 g/L NaCl as zero point.

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[] [600]					
Sample Volume Norm.	[2] [0] [0]					
Sample Volume Dec.	[20] [10] [180]					
Sample Volume Inc.	[4] [0] [0]					
Configuration:						
R1	[200] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[50] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/L]					
Decimal Places	[2]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[10]	[5]	[10]	[10]	[20]	[2]
Diluted S. Volume	[0]	[5]	[5]	[10]	[10]	[0]
Diluent Volume	[0]	[195]	[190]	[190]	[180]	[0]

* Entered by operator

3. Ordering information

AMI/AUS-000 1 x 5 mL Antiserum
 2 x 25 mL Buffer
 AMI/STD-001 AMI Standard, 1 mL
 AMI/CON-001 AMI Control, 1 mL

APPLICATION NOTE

COBAS 6000

C1 Esterase Inhibitor AUS KIT

1. Reagent preparation

Sample: Ready for use.
 Reagent 1: Buffer, ready for use.
 Reagent 3: Antiserum, ready for use.
 Calibrator: Ready for use.

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[700] [340]					
Sample Volume Norm.	[3] [0] [0]					
Sample Volume Dec.	[20] [15] [180]					
Sample Volume Inc.	[6] [0] [0]					
Configuration:						
R1	[250] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[30] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/dL]					
Decimal Places	[1]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[3]	[3]	[3]	[3]	[3]	[3]
Diluted S. Volume	[0]	[0]	[0]	[0]	[0]	[0]
Diluent Volume	[0]	[0]	[0]	[0]	[0]	[0]

* Entered by operator

3. Ordering information

AMG/AUT-000 1 x 10 mL Antiserum
 5 x 25 mL Buffer
 MPS/STS-5X1 Protein Standard Set, 5 x 1 mL
 MPS/STH001 Protein Standard High, 1 mL
 139F003 Immunology Control Low, 1mL
 139F002 Immunology Control High, 1mL

APPLICATION NOTE

COBAS 6000

Cystatin C KIT

1. Reagent preparation

Sample: Ready for use.
 Reagent 1: Buffer, ready for use.
 Reagent 3: Latex, ready for use.
 Calibrator: Cystatin C Standard, ready for use. Use 9 g/L NaCl as zero point.

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[] [600]					
Sample Volume Norm.	[2] [0] [0]					
Sample Volume Dec.	[20] [10] [180]					
Sample Volume Inc.	[4] [0] [0]					
Configuration:						
R1	[230] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[40] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/L]					
Decimal Places	[2]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[2]	[5]	[5]	[10]	[20]	[2]
Diluted S. Volume	[0]	[5]	[10]	[10]	[10]	[0]
Diluent Volume	[0]	[195]	[195]	[180]	[180]	[0]

* Entered by operator

3. Ordering information

143C001 1 x 8 mL Latex
 2 x 25 mL Buffer
 143E001 Cystatin C Standard, 1 mL
 143F003 Cystatin C Control Low, 1 mL
 143F002 Cystatin C Control High, 1 mL

APPLICATION NOTE
COBAS 6000
Kappa Light Chain AUS KIT
Serum

1. Reagent preparation

Sample: Ready for use.
 Reagent 1: Buffer, ready for use.
 Reagent 3: Antiserum, ready for use.
 Calibrator: Protein Standard High, ready for use. Use NaCl 9g/L as zero point. Alternatively, use the ready for use Protein Standard Set.

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[700] [340]					
Sample Volume Norm.	[2] [0] [0]					
Sample Volume Dec.	[20] [10] [180]					
Sample Volume Inc.	[4] [0] [0]					
Configuration:						
R1	[350] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[75] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/dL]					
Decimal Places	[1]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[2]	[5]	[10]	[20]	[20]	[2]
Diluted S. Volume	[0]	[5]	[5]	[5]	[10]	[0]
Diluent Volume	[0]	[195]	[190]	[180]	[180]	[0]

* Entered by operator

3. Ordering information

KAP/AUS-000 1 x 5 mL Antiserum
 2 x 25 mL Buffer
 MPS/STH-001 Protein Standard Set, 1 mL
 MPS/STS-5X1 Protein Standard Set, 5 x 1 mL
 139F003 Immunology Control Low, 1mL
 139F002 Immunology Control High, 1mL

APPLICATION NOTE
COBAS 6000
Lambda Light Chain AUS KIT
Serum

1. Reagent preparation

Sample: Ready for use.
 Reagent 1: Buffer, ready for use.
 Reagent 3: Antiserum, ready for use.
 Calibrator: Protein Standard High, ready for use. Use NaCl 9g/L as zero point. Alternatively, use the ready for use Protein Standard Set.

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[700] [340]					
Sample Volume Norm.	[2] [0] [0]					
Sample Volume Dec.	[20] [10] [180]					
Sample Volume Inc.	[4] [0] [0]					
Configuration:						
R1	[350] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[75] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/dL]					
Decimal Places	[1]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[2]	[5]	[10]	[20]	[20]	[2]
Diluted S. Volume	[0]	[5]	[5]	[5]	[10]	[0]
Diluent Volume	[0]	[195]	[190]	[180]	[180]	[0]

* Entered by operator

3. Ordering information

LAM/AUS-000 1 x 5 mL Antiserum
 2 x 25 mL Buffer
 MPS/STH-001 Protein Standard Set, 1 mL
 MPS/STS-5X1 Protein Standard Set, 5 x 1 mL
 139F003 Immunology Control Low, 1mL
 139F002 Immunology Control High, 1mL

APPLICATION NOTE
 COBAS 6000
Microalbumin KIT 2nd Generation
 Urine

1. Reagent preparation

Sample: Centrifuged urine
 Reagent 1: Buffer, ready for use
 Reagent 3: Antiserum, ready for use
 Calibrator: Ready for use

2. Instrument setting

Analyse						
Assay/ Time / Point	[Endpoint2] [10] [30,70,0,0]					
Wavelength (Sec./Pri.)	[700] [340]					
Sample Volume Norm.	[12] [0] [0]					
Sample Volume Dec.	[6] [0] [0]					
Sample Volume Inc.	[24] [0] [0]					
Configuration:						
R1	[170] [0] [Inactive]					
R2	[0] [0] [Inactive]					
R3	[30] [0] [Inactive]					
Dilution Water						
Dilution Diluent	[X]					
Linearity Limit	[0] [0] [0] [0]					
Prozone Limit	[-32000] [32000]					
Abs Limit	[32000] [Higher]					
Cell Detergent	[Detergent 1]					
Stirring Level	[2]					
Stirring Setting						
Calibration						
Calibration Type	[RCM]					
Point	[6]					
Span	[4]					
SD Limit	[999]					
Duplicate Limit	[20]% [100]D.O.					
Sensitivity Limit	[-99999] [99999]					
S1 Abs. Limit	[-32000] [32000]					
Range						
Unit	[mg/L]					
Decimal Places	[1]					
Automatic Rerun	[X]					
Normal Values	[*]					
Other						
Code Standard	[*]					
Concentration	[*]					
Sample Volume	[12]	[10]	[20]	[3]	[6]	[12]
Diluted S. Volume	[0]	[3]	[3]	[0]	[0]	[0]
Diluent Volume	[0]	[190]	[180]	[0]	[0]	[0]

* Entered by operator

3. Ordering information

102C002 1 x 10 mL Antiserum
 5 x 25 mL Buffer
 MAL/STD-001 Microalbumin Standard, 1 mL
 MAL/CON-001 Microalbumin Control, 1 mL
 102F003 Microalbumin Control Low, 1 mL