

Quality Control

Certificate of Analysis

QC release date: 2019-09-17

Product: Total Protein UC FS
Product Code: 10210

Cat. No.: 1 0210 99 10 021

Delivery of kits: No of kits: 25
Kit Lot/Los: 60134403
Kit Expiry: 2021-08

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 27491
Expiry: 2021-09

Examination:

Visual check: Reagent : clear, orange-red
Smell Reagent : no characteristic odour
pH: Reagent (25 °C ± 1 °C): 2.57
Absorbance: Reagent : 0.083
600 nm/ 15-25°C

Control sera: TruLab Urine 1 TruLab Urine 2
Lot: 25444 Lot: 25445
Fresh (at 37°C) 0.080 g/L 0.618 g/L
Target values 0.082 0.063 - 0.102 0.644 0.490 - 0.799

Additional remarks: Mf-Date: 03.09.2019

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka