

Quality Control

Certificate of Analysis

QC release date: 2020-01-17

Product: Lp(a) 21 FS
Product Code: 17139

Cat. No.: 1 7139 99 10 930

Delivery of kits: No of kits: 1
Kit Lot/Los: 60139353
Kit Expiry: 2021-12

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 28239
Expiry: 2021-12-31

Reagent 2 Lot No.: 28239
Expiry: 2021-12-31

Examination:

Appearance: Reagent 1: clear to slightly opalescent,
yellowish
Reagent 2: milky, white

Odour: Reagent 1: no characteristic odour
Reagent 2: no characteristic odour

Control sera: TruLab Lp(a) 1 TruLab Lp(a) 2
Lot: 27295 Lot: 27296

Test results 43.8 nmol/L 119 nmol/L

Target values 40.7 32.6 - 48.8 116 92.8 - 139

Additional remarks: Mf-Date: 02.12.2019

This lot meets our specifications.

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka